




Mirage Micro™
NASAL MASK

Mirage Activa™ LT
NASAL MASK


Mirage™ SoftGel
NASAL MASK




SMALL




MEDIUM



LARGE



LARGE WIDE



X-LARGE
(Micro only)



Actual Size: 5 1/2" (Width) x 7" (Height)

Ruler is for scaling purposes.

How to use the fitting template:

1. Cut out the template along the dotted line
2. Cut out nasal holes
(hint: fold each nasal hole in half to cut)
3. Follow directions on template to determine appropriate cushion size